Neo Line

Winkelstücke/Contra-angles/Contraᄁᄁoli
Contre-angles/Contraángulos

SL411E

Bedienungsanleitung - Operation manual - Istruzioni per l'uso
Instructions pour l'utilisation-Instrucciones de uso

Deutsch: seite 1
English: page 24
Italiano: pagina 46
Français: page 69
Español: página 92

Deutsch

Inhalt

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English

Contents

Page	Chapter
26	Information for the user
26	Symbols
27	Identification of warnings
28	Safety precautions
31	Product description
32	Purpose
33	Technical data
34	Transport and storage conditions
34	Use
39	Hygiene and maintenance
44	Disposal
44	Technical Assistance
44	Certification
45	Warranty

Information for the user

Dear Customer,

Thank you for purchasing this product from Ivory. Your new product is designed and made to the highest standards to ensure excellent performance and easy operation and maintenance. Please read the operation manual thoroughly before using the product, and keep the manual where it is readily available to all users of the product.

This document is intended for medical personnel, dentists and dental assistants.

Symbols

	Danger symbol.
	Warning symbol.
	Caution symbol.
	Hint symbol.
	Read usage instructions symbol.
	CE Mark (European Community). Products bearing this mark conform to applicable EU Directives.

	Class II medical device
	Applied Part: Type B
	Manufacturer symbol.
	This product is Autoclavable.
	This product can be washed via Thermo Disinfector.
	Do not dispose of with domestic waste.
	Request for assistance

Identification of warnings

Pay particular attention to the following warnings in order to avoid injury to persons and damage to equipment.

DANGER!

Indicates a risk that can cause death or serious injury if not avoided.

WARNING!

Indicates a risk that can cause serious injury or moderate injury if not avoided.

CAUTION!

Indicates a risk that can cause damage to equipment or moderate or slight injury if not avoided.

HINT!

General explanations, without risk to persons and/or objects.

Safety precautions

WARNING!

In the event of lesions to soft tissue in the oral cavity, or recent tooth extractions, do not perform treatments using instruments powered by compressed air. Infected substances could enter the lesions and cause infection. There is also a risk of air embolism.

CAUTION!

After treatment, remove the bur or stone before placing the handpiece back in its support. Take care not to drop!

CAUTION!

Electric micromotors produce more power than air turbines and conventional air driven motors. When handpieces are subjected to high torques and speeds, if they have received inadequate maintenance, are damaged or are used for

purposes other than that for which they are designed, they can overheat and cause serious burns to the patient. Stop using the handpiece immediately if you notice unusual functioning, excessive vibrations, overheating or any other indications of a malfunction in the device. Under normal working conditions, the surface of the handpiece must not exceed 20°C from the ambient temperature, as required by the IEC 80601-2-60:2012 standard.

WARNING!

Do not use this medical device as a light probe. Rotating burs or stones can cause oral lesions.

CAUTION!

Clean, service and dry the instrument as instructed before extended periods of disuse.

! CAUTION!

Always perform regular functioning checks and maintenance even after the warranty has expired. The following personnel are authorised to undertake the repair and maintenance of this Ivory product:

- Technicians working in Ivory factories around the world.
 - Technicians who have completed Ivory technical training.
- Technical assistance must only be performed by qualified Ivory technicians, who can also provide original Ivory spare parts.

! CAUTION!

Do not attempt to disassemble the product or tamper with any mechanism. This medical product must be used in the ways illustrated in this manual. The products and systems specified in the manual for cleaning and maintenance must also be used. Users are responsible for regular functioning checks, maintenance and inspection of the product.

! WARNING!

Sterilise the handpiece and all accessories before use.

! CAUTION!

Dirty or humid cooling air can cause malfunctioning and rapid bearing wear. Cooling air must be clean, dry and uncontaminated, and must conform to ISO 7494-2 standard.

! CAUTION!

Only use clean, undamaged burs to prevent dirt from penetrating into the locking system.

⚡ DANGER!

A bur that is loose, unsuitable or improperly seated can work free of the head or break. Before use always confirm that the bur is correctly seated and locked into place.

! CAUTION!

Insufficient cooling can lead to overheating of the

preparation and even damage the tooth. Make sure that the spray water is at least 50ml/min at a pressure of 2 bar and that the spray air is at least 1.5Nl/min at a pressure of 2 bar.

 **WARNING!**

Do not use the handpiece to move the patient's cheek while the device is functioning!

 **WARNING!**

Do not look at the LED light for a long time when it is illuminated.

 **WARNING!**

Always check the handpiece for damage before using it.

 **WARNING!**

Make sure that the micromotor in which the handpiece will be inserted has a maximum rotation of 40,000 rpm.

Before use, make sure that all the o-rings are present on the micromotor coupling.

 **WARNING!**

It can be dangerous to use the handpiece without a speed and direction controller.

- The dental unit with which this handpiece is used must be equipped with a speed and direction controller.
- For reasons of legal responsibility for safety, reliability and performance, a statement to this effect must also appear in the documentation provided with the dental unit.

 **DANGER!**

Always wear gloves to avoid the risk of lesions and infections when checking, cleaning, repairing or extracting teeth.

 **CAUTION!**

Use of unfiltered, hard water can lead to blockages and premature wear in hoses, fittings and spray heads.

Product description

Silk Line low speed contra-angle handpiece, push button, internal spray, no optic fibre, ratio 1:1. **REF. SL411E**

Purpose

This medical device is designed exclusively for the treatment of dental lesions and pathologies in a dental clinic, and conforms with the corresponding provisions of the EN ISO 14457 standard. All other uses and modifications to the product are strictly forbidden, and can lead to hazardous situations.

This medical product is intended for the following applications:

- Models SL411 E: prophylactic and restorative treatments.

All uses of this device other than those intended, and the use of accessories and spare parts not authorised by Ivory is strictly forbidden and can even be dangerous.

Failure to conform to the instructions provided herein invalidates the warranty and exonerates the manufacturer from any other claims. This medical device is to be used

only by dentists, licensed dentistry surgeons in compliance with the following:

- provisions governing safety in the workplace;
- provisions governing accident prevention;
- the instructions provided in this manual.

The user must:

- use only perfectly functional equipment;
- use the product only for its intended purpose;
- protect himself, patients and third parties against potential hazards;
- avoiding contamination by the product.

Technical data

Product code	SL411E
Indications	Prophylaxis
Ratio	1:1
Light intensity up to (lux)	/
Locking mechanism	Button
Type of spray	Internal
Type of coupling	Intramatic ISO 3964
Rotation speed (rpm) (micromotor speed set to 40,000 rpm)	40,000
Torque (Ncm)	/
Removable head	Yes

Rotatable head	No
Head diameter (mm)	9.5
Head height (mm)	11.6
Handpiece weight (g)	79
Bur/file diameter (mm)	2.334/2.350
Bur/file length (mm)	9-34
Type of lubrication	Oil
Autoclave resistant	Yes, max 135°C
Packing size lxdxh (mm)	160x90x43
Packing weight appr. (g)	250
Certification	CE 0476

Transport and storage conditions

	Temperature: from -10 °C to +55 °C (from 50 °F to +131 °F).
	Relative humidity: from 10 % to 90 % non-condensing.
	Atmospheric pressure: from 700 hPa to 1060 hPa (from 10 to 15 psi).
	Protect against humidity.
	Fragile.

Use

CAUTION!

Read the safety precautions before starting to use the instrument.

CAUTION!

At the start of every working day and after the treatment of every patient, the water-conducting systems should be rinsed for at least 2 min. without the instrument being attached.

CAUTION!

Do not disconnect the handpiece from the micromotor while it is functioning.


DANGER!

Always wear gloves to avoid the risk of lesions and infections when checking, cleaning, repairing or extracting teeth.

Fitting and removing the head


- 1) Turn the locking ring (1) fully anti-clockwise and push in the head (2) as far as it will go.
- 2) Release the locking ring (1).
- 3) Check that the head is securely held in the handpiece.
- 4) Turn the locking ring (1) fully anti-clockwise to remove the head (2).

Fitting and removing the handpiece

ONLY use this contra-angle handpiece with type E micromotors (ISO 3964).

- 1) Push the handpiece on to the micromotor until it clicks into place. In the case of handpieces with optic fibre, turn the handpiece until you hear a click, indicating correct alignment.
- 2) Check that the handpiece is securely connected to the micromotor.
- 3) Wait for the micromotor to stop before disconnecting the handpiece from it.

Fitting and removing the rotary bur

CAUTION!

Never use burs that are bent, deformed, worn, corroded, broken or defective, burs with cracks in the edges or shafts, burs that do not conform to ISO standards, burs that have been tampered with, or burs that do not conform to the specifications given herein. Follow the instructions provided by the manufacturer of the bur you are using. Only use burs that satisfy the following criteria:

- Bur shaft diameter: 2.334/2.350 mm, type 1 according to ISO 1797:2017, max. length 34 mm, code 4 according to ISO 6360-1.

The turbines are not part of the device. It is recommended to use turbines conforming to the ISO 10993 standard.

DANGER!

Always wear gloves to avoid the risk of lesions and infections when checking, cleaning, repairing or extracting teeth.

DANGER!

Only fit the bur in the handpiece when the handpiece is stopped. Never attempt to grip a bur while it is rotating, even when it is nearly stopped.

DANGER!

Never press the locking button on the handpiece while the instrument is in use or has not quite stopped. Doing so can cause the bur to become loose or cause the locking system to overheat and break.

DANGER!

AVOID CONTACT between the patient's cheek and the locking button during use. This could lead to the bur being accidentally released, resulting in burns or lesions to the patient.

WARNING!

Never operate the handpiece without a bur in the spindle.

This could cause the handpiece to overheat, causing damage or burns. Never exceed the maximum speed specified by the bur manufacturer. Keep bur shafts clean at all times; dirty shafts can cause dirt to enter the spindle, leading to the bur slipping during rotation, or even preventing the bur from fitting in the spindle correctly. Never attempt to fit an instrument to a rotating micromotor.


⚠ DANGER!

Never apply excessive pressure to burs. They could break, bend or become difficult to extract.

⚠ DANGER!

During use, avoid contact between the cheek tissues and the button; the button may be pressed inadvertently causing burn injury to the patient.


1) Press the locking button (1) and push the bur (2) in as far as it will go. On blue and green contra-angle handpieces, turn the bur until it locks in place.


⚠ WARNING!

When inserting a turbine, make sure it is inserted up to the stop. Before insertion, make sure the shank of the turbine is clean.

2) Check the tightness of the turbine (2).

3) Press the locking button (1) to release the bur (2).

⚡ DANGER!

If you notice any unusual noise, vibration, overheating, interruption in the flow of cooling liquid, leakage, stop using

! CAUTION!

Always check that the cooling system is unobstructed before starting work.

Hygiene and maintenance

⚠ WARNING!

The maintenance processes are described according to ISO 17664. Before starting a process:

⚡ DANGER!

Always wear gloves to avoid the risk of lesions and infections when checking, cleaning, repairing or extracting teeth.

! CAUTION!

Never immerse the handpiece in disinfectant solutions or ultrasonic baths! Only use chlorine-free disinfectant products that have been certified by officially recognised institutes.

- Manual external cleaning

Clean and disinfect the handpiece after every treatment to eliminate residues of bodily fluids (e.g. blood, saliva etc.) and

prevent them entering the internal parts of the handpiece.


! CAUTION!

The saline solution is extremely aggressive and can rust steels in no time.

Thoroughly clean the surfaces using a clean cloth soaked in a suitable substance. Also use a medium hard toothbrush and brush under running hot tap water $30^{\circ}\text{C} \pm 5^{\circ}\text{C}$ ($86^{\circ}\text{F} \pm 10^{\circ}\text{F}$).

- Cleaning the spray nozzles

The use of hard, unfiltered water prematurely clogs the tubes and spray diffusers. To remove any encrustations in the three spray nozzles, use the special probe. Using the supplied cleaning wire, carefully remove impurities and deposits from the spray nozzles. Clean and disinfect the cleaning thread before storing it.


! CAUTION!

Do not force the brush into the conduit. Damage to the conduit can change the direction of spray with respect to the bur and reduce the effectiveness of cooling.

- Cleaning the optic fibre glass (if fitted)

Clean the glass of the optic fibre with a pad dampened in alcohol. Remove all residual impurities.

Manual internal cleaning

Using a suitable spray product, follow the corresponding instructions and spray the substance into the medical device. Immediately after internal cleaning, carry out the disinfection process or the drying process.

- Automatic cleaning

Use an ISO 15883 compliant washer-disinfector that uses an alkaline cleaning agent with a maximum pH of up to 10.

- Disinfection

External manual: disinfect using a clean cloth soaked in a suitable substance (e.g. O-phenylphenol or ethyl alcohol). Products not recommended: disinfectants based on benzalkonium chloride, acetone or glutaraldehyde.

Internal manual: Using a suitable spray product, follow the corresponding instructions and spray the substance inside the medical device. Immediately after internal disinfection, perform the oil lubrication process and / or the sterilisation process.

Automated: Use an ISO 15883 compliant washer-

disinfector that uses an alkaline cleaning agent with a maximum pH of up to 10.

- Drying

Manual: in general use dry and clean air compliant with ISO 7494-2.

Automated: the drying phase is normally part of the cleaning programme of the washer-disinfector.

- Lubrication

 **CAUTION!**

Lubrication cycles:


- Before every sterilisation.
- After every 15 minutes of rotation.

 **CAUTION!**


Lubricate only with Ivory oil OL101S.

Procedure to be performed:


1) Turn the locking ring (1) fully anti-clockwise and remove the head.

2) Fit the adapter for Ivory heads and contra-angles to the oil spray can. Fit the head (A) in the adapter (B). Hold the head in place. Spray for about 5 seconds.

3) Apply the adapter to the back of the handpiece. Hold the handpiece in place and spray for about 5 seconds.

4) Replace the head in the handpiece.

! CAUTION!

Incorrect or inadequate spindle cleaning can cause the spindle to loosen and accidentally release the bur during functioning.

! CAUTION!

After lubricating, hold the handpiece with the chuck facing down. Operate the handpiece for 30 seconds to expel any excess oil. If you notice any dirt coming out, repeat the entire cleaning and lubrication process.

-Sterilisation

Sterilisation in steam autoclave with fractional pre-vacuum, class B cycle according to the standard: EN 13060, EN ISO 17665-1

CAUTION!

- Follow the instructions of the autoclave manufacturer.
- Clean, disinfect and lubricate the handpiece before sterilisation. The presence of blood residues on the external or internal surfaces can lead to the formation of corrosion.
- Make sure that the handpiece is properly dried after sterilisation.

Authorised sterilisation procedures

- 1) Seal the handpiece and accessories in bags compliant with UNI EN 868:2009, UNI EN ISO 11607:2009 "Packaging for terminally sterilised medical devices".
- 2) Sterilise the bagged handpiece in an autoclave according to **EN ISO 17665-1** standard and respect the following conditions:

- Class B steam sterilisation in autoclaves conforming to **EN 13060** standard. Sterilisation time at least 3.5 minutes at 135 °C.
- 3) Remove the bag from the autoclave.

-Storage

- Keep the handpiece in the sterilisation bag until it is needed.
- Once sterilised, store the product in a dry and dust-free place, if possible with a low microbial charge. The processed device must be stored protected from dust and microbiotic contamination. Storage temperature -10 + 70°C. Storage humidity 10 + 90%.

ATTENTION!

Observe the validity date of the sterilisation carried out!
The device is supplied "non-sterile". The device can be sterilised for at least 1000 times.

Disposal

 Sterilise the product and its accessories before disposing of them. Refer the disposal of the product to authorised companies qualified to dispose of special industrial waste.

Technical assistance – Ivory direct line

 **HINT!**

To request technical assistance, compile the on-line form available at: www.ivory-dent.com

Declaration of conformity

This medical device conforms to the requirements of European Directives 93/42 CEE (Class IIa) and 2007/47/CE. This conformity is expressed by the mark:

CE
0476

Additional documentation for controller n°0476 can be requested from: ivory@ivory-dent.com

Warranty

Ivory guarantees this product against defects of material and manufacture for a period of 24 months from the date of purchase. Ivory cannot be held responsible for defects and consequences caused by or attributable to normal wear, improper use, cleaning or maintenance, failure to follow operating and coupling instructions, lime scale or corrosion, contamination of the air or water supplies, chemical or electrical effects of an exceptional nature or not permitted by the instructions provided by Ivory or other manufacturers. Ivory expressly declines all responsibility for defects or consequences attributable to tampering with the product or modifications made by the customer or unauthorised third parties, or attributable to the use of non-original spare parts. In the event of justified claims, Ivory shall either repair the product or supply spare parts free of charge under the warranty. The warranty excludes all other claims, with particular reference to claims for compensation. In the case of delays, serious or wilful negligence, this exclusion is only valid provided no specific laws are contravened.

The warranty does not apply to lamps, glass and optic fibre conductors, glass or rubber parts, or the colour fastness of plastic materials or all parts subject to normal wear. Claims made under the warranty must be presented by compiling the on-line request "QRRMA" form and by accompanying the product with a copy of the purchase invoice.

The purchase invoice must clearly identify the dealer, date of purchase, code or model of the product and the production or serial number. Please refer also to the articles 17,18 and 19 in the "Terms and payment" section on www.ivory-dent.com

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KaVo® is a registered trademark of Kaltenbach & Voigt GmbH & Co., Germany. Sirona® is a registered trademark of Sirona Dental Systems GmbH, Germany. Bien-Air® is a registered trademark of Bien-Air Dental S.A., Switzerland. NSK® is a registered trademark of NSK NAKANISHI INC., Japan. W&H® and Roto Quick® are registered trademarks of W&H Dentalwerk Bürmoos GmbH, Austria.